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1. GENERAL INFORMATION

The information on this information sheet refers to the use of emerging monolithic implants Ø 2.8 - 3.5 - 4.1 - 5.0 mm, with different lengths.

All Bone System® products are designed and manufactured as an essential part of a structured, composite implant system. When fitting the system, use original Bone System® components, accessories and instruments only, carefully following the manufacturer's instructions and advice.

In accordance with applicable regulations, all Bone System® products undergo stringent quality control tests, to provide a flawless product that is free of obvious defects.

1.1 Warnings

- Read the instructions for use carefully and familiarise with the components and instruments before each use.
- The instructions provided on this sheet cannot guarantee perfect results when the products are used by practitioners who are not specialised in implantology. Consequently, both the implant surgeon and the prosthetic technician must be experts in oral surgery and dental prosthetics, including diagnosis, planning and lab procedures. It is therefore highly recommended that both the implant surgeon and the prosthetic technician take part in the theoretical and practical training courses that Bone System® organises for its customers.
- Bone System® products are intended for use only by dentists and dental technicians, working within their area of competence and in compliance with safety regulations.
- All the products described herein should be used in compliance with the principles of dentistry. Incorrect assessment of the case, from both a clinical and a biomechanical and prosthetic point of view, could result in treatment failure.
- Users apply, machine and handle each product according to their autonomous judgement and consequently assume full responsibility for their work. No responsibility may be attributed to Bone System® for damage of any kind caused by such actions.
- If you have any doubts, do not use the products without first seeking the necessary clarifications.

1.2 Precautions for use

- Before performing any procedure, make sure you are in possession of the necessary amounts of all the components and appropriate instruments in good working order.
- Use original components and instruments in perfect condition only.
- Position the patient so as to avoid the risk of him/her swallowing or inhaling the components during use.
- Prepare the implant site following the protocol provided for each type of implant.
- Always remember that the tip of the drill reaches a greater depth than that required for the implant. This is very important in close proximity to anatomical elements at risk that must be preserved and in such cases an adequate safety margin (1-2 mm) must be respected. An error in the evaluation of the depth during implant site preparation could cause permanent nerve damage or cause paraesthesia or other, no less serious, injuries.
- **During the preparation of the implant site and the insertion of the implants, they must be positioned parallel to natural teeth and adjacent structures.**
- When drilling inside the mouth, use suitable irrigation and aspiration.

1.3 Description and instructions for use

Bone System® emerging monolithic implants "QM" incorporate into a single component the endosseous portion, the transmucosal height and the abutment. This makes them particularly well-suited to minimally-invasive and flapless techniques, when appropriate to the patient's clinical situation. These implants can be used in restoration work on both arches, in all cases of single, uni- or bilateral tooth loss and in completely edentulous patients. They are appropriate for immediate loading techniques and to support fixed temporary or permanent cement-retained prosthetic structure.

WARNING! - The use of the Ø 2.8 mm implant is limited to the front areas. The Ø 2.8 mm implant is ideal for mono-edentulous procedures with limited mesial-distal space or in combination with other implants required by the treatment plan.

1.4 Contraindications and side effects

Implantology articles and manuals provide an exhaustive overview of patient assessment criteria, the absolute or relative contraindications of implant treatment and the most suitable surgical techniques for each case.

The main absolute contraindications for implant surgery are:

- patients who have received radiation therapy in the past year and those who are receiving cancer treatments or treatments containing bisphosphonates;
 - uncontrolled diabetes;
 - over SLE;
 - mental imbalance;
 - genetic factors making the subject prone to bone loss (genetic polymorphism).
- When used correctly, there are no known specific contraindications associated with the use of Bone System® implant system components.

1.5 Materials

Bone System® emerging monolithic implants "QM" are manufactured in titanium, with physical and mechanical characteristics suited to their use, in accordance with applicable regulations.

1.6 Packaging

Bone System® emerging monolithic implant system "QM" are packed in blister packs and sterilized by gamma rays in accordance with European regulations. Successful sterilization is proven by a ball sticker on the pack that turns red after the treatment.

If preserved, the blister pack isolates the content from external agents and guarantees its sterility (up to the expiry date) until the pack is opened, unless the blister is damaged.

1.7 Storage

Bone System® emerging implant system "QM" must be stored in their original packs and in a dry place such as to maintain their condition.

2.1 Opening the pack

The sterile pack should be opened immediately before the introduction of the implant and when handling the content, the normal rules of asepsis must be applied to avoid contamination.

The implant pack bears a label providing the information needed for identification. A similar label, inside the pack, must be attached to the patient's medical records for implant traceability purposes.

2. MODE OF USE


2.2 Picking up and inserting the implant

1. Open the blister and turn the tube out onto the sterile field

2. Remove the cap from the tube together with the connected implant and introduce into the patient's mouth, placing the implant inside the previously prepared site.

3. Using the cap attached to the implant, slowly and gradually screw the implant into the site, until you feel a certain resistance.

4. When the implant is inserted and adequately stable, remove the cap and complete the screwing operation with a specific driver, fitted on a manual knob tool or a dynamometric ratchet. The implant's special characteristics mean that if necessary, it is possible to adjust the angle at which it is inserted in relation to the direction of the implant site.

NOTE - Given the geometry of the implant, and relative preparation of the implant site, the implant must be inserted into the bone so that only the nitride portion (transmucosal height + abutment) protrudes.

WARNING! If difficulties are encountered during the insertion of the implant, loosen the implant a few turns and try again. If the problem persists, reconsider the adequacy of site preparation. Do not exceed torque values of 45-50 Ncm when tightening the implant to avoid damaging the bone and compromising the conditions of the tools used.

2.3 Parallel placement and preparation of the emerging abutment


IMPORTANT NOTE - Bone System® emerging monolithic "QM" implants have been designed to rectify slight divergences inside the patient's mouth by drilling or folding of the protruding part. The parallelism procedure is easiest to perform using Ø 2.8 and 3.1 implants. It is important to remember that the risk of fracturing the cortical bone in the implant site increases with an increase in implant diameter.

- Perform folding before suturing the flaps (unless a flapless technique has been chosen) to avoid damaging the soft tissues around the fold area.

- Using the parallelism tool, work the protruding abutment very carefully to give it the required direction. During this step, grasp the alveolar crest between the fingertips to prevent any cortical bone collapse.

- Proceed with very small, gradual movements in the same direction, avoiding movements that could warp and weaken the abutment, with a consequent risk of fracture in the future.

- After suturing, isolate the abutment with a rubber dam; adjust the height of the abutment and fine-tune the protruding portion, if necessary. These operations should be performed using a bur suitable for metal, under abundant irrigation with chilled saline solution, aspirating the metal residues produced by the drilling

3. ADDITIONAL INFORMATION

Bone System® products are identified according to the instructions and symbols pursuant to European Directive EN 15223, namely:

1. DO NOT reuse. DISPOSABLE device, to be used only once for one patient. Reuse could involve serious risks of cross-infection and/or functional loss.
2. Use before (expiry date)
3. Batch code
4. Sterile - sterilized by ionising radiation

5. NOTSterile
6. DO NOT use if the pack is not intact
7. Product code
8. Warning, see supporting documents
9. Manufacturer
10. See the instructions for use

Not all of the symbols mentioned are present on the pack of the product.

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All mentioned trademarks are the property of their respective owners. © Bone System. The copying of all or part of the text and/or images present on this information sheet is prohibited unless explicitly authorized.