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

1. GENERAL INFORMATION

These instructions concern the following prosthetic and auxiliary components:

- A) Impression abutments
- B) Temporary screw-retained abutment
- C) Abutment for the Easy-on-Four® technique®
- D) Straight or offset screw-retained Friction Locking® abutment
- E) Abutment for the Tissue level® technique
- F) Impression transfer and screw for the Easy-on-Four® technique
- G) Prosthetic cap and fixing screw for the Easy-on-Four® technique

All Bone System® products are designed and manufactured as an essential part of a structured, composite implant system. When fitting the system, use original Bone System® components, accessories and instruments only, carefully following the manufacturer's instructions and advice.

In accordance with applicable regulations, all Bone System® products undergo stringent quality control tests, to provide a flawless product that is free of obvious defects.

1.1 Warnings

- Read the instructions for use carefully and familiarise with the components and instruments before each use.
- The instructions provided on this sheet cannot guarantee perfect results when the products are used by practitioners who are not specialised in implantology. Consequently, both the implant

surgeon and the prosthetic technician must be experts in oral surgery and dental prosthetics, including diagnosis, planning and lab procedures. It is therefore highly recommended that both the implant surgeon and the prosthetic technician take part in the theoretical and practical training courses that Bone System® organises for its customers.

- Bone System® products are intended for use only by dentists and dental technicians, working within their area of competence and in compliance with safety regulations.
- All the products described herein should be used in compliance with the principles of dentistry. Incorrect assessment of the case, from both a clinical and a biomechanical and prosthetic point of view, could result in treatment failure.
- Users apply, machine and handle each product according to their autonomous judgement and consequently assume full responsibility for their work. No responsibility may be attributed to Bone System® for damage of any kind caused by such actions.
- If you have any doubts, do not use the products without first seeking the necessary clarifications.

1.2 Precautions for use

- Before performing any procedure, make sure you are in possession of the necessary amounts of all the components and appropriate instruments in good working order.
- Use original components and instruments in perfect condition only.
- Position the patient so as to avoid the risk of him/her swallowing or inhaling the components during use.
- When drilling inside the mouth, use suitable irrigation and aspiration.
- Tighten the screws to the values indicated to avoid excessive torque or loosening due to inadequate tightening.

1.3 Intended use

These components are compatible for use with all Bone System® implant system with internal connections, designed for impression-taking and prosthetic restoration work, as described in section 2.

1.4 Contraindications and side effects

There were no known contraindications at the time this instruction sheet was compiled.

1.5 Materials

All components mentioned in this instruction sheet are manufactured in titanium, with physical and mechanical characteristics suited to their use.

1.6 Packaging

The prosthetic components described here are supplied NON-STERILE in a medical bag NOT SUITABLE for autoclaving.

1.7 Storage

Products must be stored in their original packs and in a dry place such as to maintain their condition.

2. MODE OF USE

WARNING: NON-sterile components: Wash and sterilize before use according to its own previously validated method.

2.1 Impression abutments (A)

- Select the abutment best suited to the height of the collar or collarless situation.
- Insert the hexagonal part of the abutment inside the implant, pushing it all the way in.

2.2 Temporary screw-retained abutment (B)

- Select the abutment best suited to the height of the collar or collarless situation.
- Tighten by hand using a 2.3 square driver.
- Shorten or modify the protruding portion as necessary.

2.3 Straight and offset abutments for the Easy-on-Four® technique (C)

- Select the straight or offset abutment that best suits the clinical situation.
- Type 1** - Tighten straight abutments with the 3.2 hexagonal driver and the dynamometric ratchet at a maximum torque of 32 Ncm.
- Type 1** - Tighten offset abutments using the 1.27 hexagonal driver and the dynamometric ratchet at a maximum torque of 32 Ncm.

2.4 Straight and offset screw-retained Friction Locking® abutments (D)

- Select the straight or offset abutment according to the clinical situation.
- Tighten the straight or offset abutments with the 1.27 hexagonal driver and dynamometric ratchet at a maximum torque of 32 Ncm.

2.5 Tissue level® abutment technique (E)

- Type 1** - Tighten the abutment with the "D-Shape" driver and the dynamometric ratchet at a maximum torque of 32 Ncm.
- Type 2** - Tighten the abutment with the 1.27 hexagonal driver and the dynamometric ratchet at a maximum torque of 32 Ncm.

2.6 Impression transfers and screws for Easy-on-Four® technique (F)

- Apply the transfer and screws to the abutment
- Tighten the transfer screws by hand until they are fully inserted.

2.7 Titanium prosthetic cap and screw for Easy-on-Four® technique (G)

- Adjust the caps to the clinical situation and make the prosthetic structure to incorporate them correctly.
- Position the prosthetic structure with incorporated caps inside the mouth.
- Tighten the fixing screws using the 1.27 hexagonal driver and dynamometric ratchet at a maximum torque of 22-25 Ncm.

3. ADDITIONAL INFORMATION

Bone System® products are identified according to the instructions and symbols pursuant to European Directive EN 15223, namely:

- | | |
|--|--|
| | 1. DO NOT reuse. DISPOSABLE device, to be used only once for one patient. Reuse could involve serious risks of cross-infection and/or functional loss. |
| | 2. Use before (expiry date) |
| | 3. Batch code |
| | 4. Sterile - sterilized by ionising radiation |
| | 5. NOT Sterile |
| | 6. DO NOT use if the pack is not intact |
| | 7. Product code |
| | 8. Warning, see supporting documents |
| | 9. Manufacturer |
| | 10. See the instructions for use |

Not all of the symbols mentioned are present on the product pack.

All mentioned trademarks are the property of their respective owners. © Bone System. The copying of all or part of the text and/or images present on this information sheet is prohibited unless explicitly authorized.

Bone System s.r.l.
I-20134 MILANO - ITALY - Via Rubattino 94/A
Tel. 02.21.54.165 - Fax 02.21.54.292
www.bonesystem.it
e-mail: info@bonesystem.it

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